
**Botron** Bright yellow disposable dissipative heel grounders are made from the best products on the market. These units are ideal for visitors and will cut the cost of using high priced standard heel grounders. They are one size fits all and provide excellent visibility. All visitors must wear two grounders. **Disclaimer.** All statements of technical information are believed to be true and are based upon tests we believe to be reliable. The proper use and application for this product must be the responsibility of the user. The statements herein shall have no force or effect.



Release paper with pressure sensitive adhesive



INSTALLATION

- Remover release paper
- Stick on bottom of heel of shoe
- Pull material tight up and tuck into shoe
- Step firmly onto pressure sensitive area and tear off excess material at serrations.

#### **OTHER PRODUCTS:**

- B7500 Ergo-One heel grounder
- B7571 Non-Marking heel grounder
- B7525 Snap-Loc heel grounder
- **B7591** Red Hi-Visibility heel grounder
- B7522 Toe grounder

Many other heel & toe grounders

# BOTRON CO. INC.

325 W. MELINDA LN PHOENIX AZ 85027 PH# 623-582-6700 FAX# 623-582-6776

## FINAL TEST

Calle del Ebano #16625 Jardines de Chapultepec Tijuana B.C. Mexico Tel. (664) 681 1130 Fax. (664) 681 1150 Tel. 01800 027-4848 www.finaltest.com.mx Bright yellow material with conductive fibers

PROPERTIES	SPECIFICATIONS
RTT:	5x10 <sup>6</sup>
Material:	Dissipative material
Tab:	15" long
Charge Decay:	<0.01Sec
Size:	15"
Color:	Bright yellow
Fastening system:	Pressure sensitive adhesive
Packing:	100 per. bag

#### ALWAYS USE TWO GROUNDERS

<b>RELATED TEST EQUIPMENT:</b>	
B82250	Dual Tester with stand
B82450	Combo tester with stand
B8525	Deluxe Dual tester & stand
B85502	Calibration unit for testers
B8506	Wall test station

### DISPOSABLE HIGH VISIBILITY HEEL GROUNDER

S750368